

There was little difference in endothelial cell loss (ECL) between the iStent inject and control groups. Results were consistent with previous reports of cataract surgery-related ECL. The mean percent change in ECD from baseline to 24 months was -13.1% (SD 12.4; 95% CI -14.4%, -11.8%) for the iStent inject group and -12.3% (SD 12.7%; 95% CI -14.8%, -9.8%) for the control group.

A similar proportion of eyes in each group (10.4% in the iStent inject group and 9.5% in the control group) experienced ECL > 30% at 24 months postoperatively.

2. Effectiveness Results

Results from the primary and secondary endpoints are shown in **Table 9**. The primary effectiveness endpoint was met, with 75.8% (288/380) in the iStent inject group and 61.9% (173/278) in the Control group achieving a clinically significant ($\geq 20\%$) reduction in medication-free diurnal IOP from baseline at 24 months. This difference between groups was statistically significant ($p=0.003$).

The secondary endpoint, a clinically significant mean change in medication-free diurnal IOP from baseline at 24-month postoperative examination, was met. The mean reduction in medication-free mean diurnal IOP from baseline to 24 months was 7.0 mmHg (SD 4.0) in the iStent inject group compared to 5.4 mmHg (SD 3.7) in the control group ($p < 0.001$).

Effectiveness Endpoint (Evaluated at 24 Months Postoperatively)	Cataract Surgery with iStent inject N = 380	Cataract Surgery Only N = 278	Difference (iStent inject vs. control)	P-value for difference
Proportion of subjects with medication-free DIOP reduction $\geq 20\%$ from baseline	75.8%	61.9%	13.9%	0.003 ¹
Medication-free mean DIOP (mmHg) change from baseline ²	-7.0	-5.4		< 0.001 ¹

Subjects without Month 24 medication-free diurnal IOP or with IOP-related SSIs, loss of light perception or hypotony (IOP < 6 mmHg) associated with clinically significant findings prior to 24 months were treated as non-responders. iStent inject subjects with stent reposition or removal prior to 24 months were treated as non-responders.

- The 24-month diurnal IOP values were subtracted from baseline diurnal IOP in all subjects, except for the non-responders described above. For the non-responders described above, the baseline diurnal IOP values were used for the 24-month diurnal IOP values (i.e., a change of 0 mmHg was used).
- One-sided Fisher's exact test with a significance level of 0.025.
- One-sided two-sample t-test with a significance level of 0.025.

Additional detail regarding the reasons patients did not achieve the primary endpoint (IOP non-responders) is shown in **Table 10**.

	Cataract Surgery with iStent inject n/N (%)	Cataract Surgery Only n/N (%)
Total Non-Responders	91 (24.2%)	41 (18.7%)
Non-Responders: 24-month unmedicated diurnal IOP reduction from baseline < 20%	56 (14.7%)	26 (22.0%)
Non-Responders for reasons other than IOP reduction ¹	36 (9.5%)	19 (16.1%)
Secondary glaucoma surgery ²	5 (1.3%)	3 (2.5%)
Other IOP-affecting secondary surgery ³	0 (0.0%)	0 (0.0%)
Stent reposition or removal	0 (0.0%)	0 (0.0%)
Loss of light perception	0 (0.0%)	0 (0.0%)
Clinically significant hypotony	0 (0.0%)	0 (0.0%)
Did not complete medication washout – Safety concerns	12 (3.2%)	4 (3.4%)
Did not complete medication washout – Instructions not provided/followed ⁴	0 (0.0%)	2 (1.7%)
Missing 24-month diurnal IOP data ⁴	19 (5.0%)	10 (8.5%)
Death	4 (1.1%)	6 (5.1%)
Investigator's decision	1 (0.3%)	0 (0.0%)
Last contact	0 (0.0%)	2 (1.7%)
Subject's decision	6 (1.6%)	2 (1.7%)

- n = number of eyes with the corresponding responses, % = n / N x 100%.
- Subjects were included in the primary category of "Non-Responders for reasons other than IOP reduction".
 - Secondary glaucoma surgeries include trabeculectomy, and laser trabeculectomy.
 - Other IOP-affecting secondary surgeries.
 - The outcomes of these subjects were imputed for the 24-month analysis. There were 2 subjects on oral medication at 23 months and both subjects underwent washout. Hence, although any subjects on oral medication at 24 months would have been considered non-responders due to the potential to confound the endpoint analysis, there were no subjects in this category.

- Summary of Supplemental Clinical Information
 - For the pivotal trial of the iStent inject, the Ocular Surface Disease Index (OSDI[®]) was self-administered by study subjects. The OSDI questionnaire contains 12 questions involving ocular symptoms, vision-related function and environmental triggers experienced by the subject during the past week, and is assessed on a scale of 0 to 100 with higher scores representing greater disability. **Table 11** summarizes the change in OSDI subscales and overall score from baseline. The mean improvements at 24 months from baseline were slightly higher in the iStent inject group compared to the control group involving ocular symptoms (16.41 vs. -10.69) and vision-related function (22.60 vs. -18.56) and similar involving environmental triggers (-7.41 vs. -7.70). The mean improvement in OSDI overall score at 24 months was also higher in the iStent inject group compared to the control group (-16.25 vs. -12.38). The questionnaire used to collect these data has not been validated, and therefore the true rates of these symptoms may differ from those presented in the **Table 11**.

Statistics	Cataract Surgery with iStent inject Total Number of Subjects = 386				Cataract Surgery Only Total Number of Subjects = 319			
	1M n (%)	6M n (%)	12M n (%)	24M n (%)	1M n (%)	6M n (%)	12M n (%)	24M n (%)
Ocular Symptoms (01, 02, 03)								
N	382	376	367	361	117	118	115	109
Mean	-11.87	-15.04	-16.93	-16.41	-6.41	-10.55	-11.53	-10.69
SD	22.39	21.23	19.96	21.13	20.53	18.45	17.16	17.74
Median	-10.0	-15.0	-15.0	-15.0	-5.0	-10.0	-10.0	-10.0
Min	-100	-100	-90.0	-100	-55.0	-60.0	-75.0	-65.0
Max	75.0	50.0	33.8	60.0	80.0	40.0	35.0	35.0
Not Reported	2	1	3	5	2	0	1	0
Vision-Related Function (04, 05, 06, 07, 08, 09)								
N	379	374	363	359	117	118	115	109
Mean	-16.07	-21.46	-22.82	-22.60	-14.08	-17.32	-20.92	-18.56
SD	29.80	27.93	28.22	27.30	29.94	27.49	27.66	28.92
Median	-12.5	-18.8	-18.8	-18.8	-6.3	-12.5	-16.7	-12.5
Min	-93.8	-100	-100	-100	-100	-100	-100	-100
Max	100.0	77.1	62.5	62.5	87.5	75.0	37.5	68.8
Not Reported	5	3	7	7	2	0	1	0
Environmental Triggers (10, 11, 12)								
N	370	367	358	353	114	116	113	106
Mean	-5.20	-7.27	-7.83	-7.41	-4.61	-7.26	-7.82	-7.70
SD	21.52	20.70	21.65	22.61	21.95	21.61	21.60	20.66
Median	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Min	-83.3	-100	-100	-100	-100	-100	-100	-75.0
Max	100.0	58.3	75.0	66.7	66.7	41.7	33.3	75.0
Not Reported	14	10	12	13	5	2	3	3
Overall Composite Score								
N	382	376	367	361	117	118	115	109
Mean	-11.87	-15.44	-16.66	-16.25	-8.48	-11.91	-13.60	-12.38
SD	20.79	19.39	19.38	19.73	20.02	18.01	17.18	18.38
Median	-10.4	-12.5	-13.3	-13.3	-6.7	-10.4	-10.7	-10.4
Min	-93.8	-93.8	-95.8	-100	-60.4	-66.7	-64.6	-43.5
Max	72.9	37.5	31.3	48.8	78.8	37.5	17.6	56.3
Not Reported	2	1	3	5	2	0	1	0

- Each sub-scale is a summation of some specific questions to the OSDI.
- Based on proportional analysis using a non-responder imputation for missing data. Subjects without Month 24 medication-free diurnal IOP or with IOP-related SSIs, loss of light perception or hypotony (IOP < 6 mmHg) associated with clinically significant findings prior to 24 months were treated as non-responders. iStent inject subjects with stent reposition or removal prior to 24 months were treated as non-responders.
 - In the iStent inject pivotal trial, at 24 months, the proportion of subjects with medication-free diurnal IOP ≤ 18 mmHg was 63.2% in the treatment group and 50.0% in the control group (difference: 13.2%; 95% CI 2.9%, 23.4%).
 - In the iStent inject pivotal trial, mean observed unmedicated IOP was higher at baseline and lower at 24 months in the iStent inject group. IOP at baseline was 24.8 (SD 3.4) mmHg in the iStent inject group and 24.5 (SD 3.1) mmHg in the control group. Unmedicated IOP at 24 months was 17.1 mmHg (SD 3.6) at 24 months in the iStent inject group and 17.8 mmHg (SD 3.5) in the control group¹.
 - Of the subjects who were responders (e.g., 24-month unmedicated mean DIOP was reduced by $\geq 20\%$ as compared with baseline in the absence of IOP-affecting surgery during the study), 84% of subjects in the iStent inject group (143/288) and 67% of subjects in the Control Group (49/73) were not using ocular hypotensive medication at 23 months.

¹Based on mean observed unmedicated IOP values from only those subjects with unmedicated IOP and without SSIs or other events (including loss of light perception or hypotony (IOP < 6 mmHg) associated with clinically significant findings).

14. LABELING

Not all symbols may be included in the labeling of this product.

Symbol	Definition	Symbol	Definition
	Catalogue/Model Number		Consult Instructions for Use
	Serial Number (for the stent)		Manufacturer
	Lot Number		Sterilized by Gamma Irradiation
	Do not reuse		For prescription use only
	Use-by date (year-month-day)		Temperature Storage Requirement
	Do not use if package is damaged		MR Conditional

15. MRI SAFETY INFORMATION



Non-clinical testing has demonstrated that the iStent inject W Trabecular Micro-Bypass System Model G2-W is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3 T or less
- Maximum spatial gradient magnetic field of 4.000 gauss/cm (40 T/m)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 4 W/kg

Under the scan conditions defined above, the iStent inject W Trabecular Micro-Bypass System Model G2-W is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends less than 15 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

16. CAUTION

Federal law restricts this device to sale by, or on the order of, a physician. Physician training by certified Glaukos personnel is required prior to use of this device. Training consists of three main parts:

- Didactic session
- Simulated implantation of iStent inject W
- Supervised iStent inject W implantation of clinical cases until implantation proficiency is demonstrated

Manufacturer:

Glaukos Corporation
 229 Avenida Fabricante
 San Clemente, CA 92672 U.S.A.
 Tel: +1 949 367 9600, Fax: +1 949 367 9984
 www.glaukos.com
 Toll-Free: 1-800-GLAUKOS (452-8567)

Patented: Patent info: www.glaukos.com/patents

Glaukos[®] and iStent inject[™] W are registered trademarks of Glaukos Corporation. Tyvek[®] is a registered trademark of DuPont USA.