INSTRUCTIONS FOR USE/PACKAGE INSERT
Glaukos Corporation iStent® Trabecular Micro-Bypass Stent System

DEVICE DESCRIPTION

The iStent Trabecular Micro-Bypass Stent Model GTS100R/L is an intraocular stent that is designed for implantation into subjects diagnosed with mild to moderate open-angle glaucoma. The iStent Trabecular Micro-Bypass Stent is manufactured from titanium (Ti6Al4V ELI) and is heparin coated (note: the heparin is from a porcine source). The stent has a single piece design, is 1.0 mm in length, 0.33 mm in height, with a snorkel length of 0.25 mm, and a nominal snorkel bore diameter of 120µm. The implant is provided to the surgeon in this pre-loaded configuration in order to allow for easy insertion into Schlemm’s canal. The inserter has been designed by Glaukos Corporation to hold the implant and to release the implant once it has been inserted within Schlemm’s canal. Two versions of the stent, a left-eye stent (Model GTS100L) and a right-eye stent (Model GTS100R) are available. The stents are identical, except that the “foot” faces the opposite direction in order to facilitate nasal placement and to optimize distal outflow through the collector channels.

The Glaukos iStent Trabecular Micro-Bypass Inserter Model GTS100i is also available in a stand-alone configuration; i.e., The Inserter Model GTS100i does not have an iStent Trabecular Micro-Bypass Stent attached to the tip when packaged in this configuration. The inserter is provided as a single-use, disposable device in order to allow for repositioning or removal of stents perioperatively or postoperatively. The inserter is able to reacquire the stent intraocularly should the surgeon prematurely release the stent.

Glaukos Corporation iStent Trabecular Micro-Bypass Stent System

Catalogue #	Description
GTS100L	Left eye iStent attached to disposable inserter
GTS100R	Right eye iStent attached to disposable inserter
GTS100i	Stand-alone inserter (no stent attached)

INDICATIONS FOR USE

The iStent Trabecular Micro-Bypass Stent Model GTS100R/L is indicated for use in conjunction with cataract surgery for the reduction of IOP in subjects with mild to moderate open-angle glaucoma currently treated with ocular hypotensive medication.

CONTRAINDICATIONS

The iStent Trabecular Micro-Bypass Stent is contraindicated under the following circumstances or conditions:

- In patients with angle closure glaucoma and other secondary glaucomas such as neovascular glaucoma and uveitic glaucoma.
- In patients with retrolubar tumor, chronic inflammatory disease, thyroid eye disease, Sturge-Weber Syndrome or any other type of condition where the trabecular meshwork, Schlemm’s Canal or collector channels at the implant site are compromised.

WARNINGS/PRECAUTIONS

1. The effectiveness of the iStent Trabecular Micro-Bypass Stent has not been established in patients with pseudophakic, secondary, uveitic, or neovascular glaucoma or discernible congenital anomalies of the anterior chamber angle.
2. The use of the iStent Trabecular Micro-Bypass Stent has not been studied as an alternative to the primary treatment of glaucomatous symptoms with medications. The effectiveness of this device has been demonstrated only in patients with mild to moderate open-angle glaucoma who are currently treated with ocular hypotensive medication and who are undergoing concurrent cataract surgery.
3. The surgeon should monitor the patient postoperatively for proper maintenance of intraocular pressure. If intraocular pressure is not adequately maintained after surgery, the surgeon should consider an appropriate medication regimen to reduce intraocular pressure.
4. The safety and effectiveness of the iStent Trabecular Micro-Bypass Stent has not been established in patients that experience significant complications during cataract surgery, including but not limited to, severe corneal burn, vitreous removal/vitrectomy required, corneal injuries, or complications requiring the placement of an anterior chamber IOL.

INSTRUCTIONS FOR USE

1. Cataract Surgery

- a. Preoperatively, the subject should be directed to use a fluoroquinolone antibiotic 4 times per day for one day preoperatively **in the operative eye**. An additional drop should be administered 30 minutes preoperatively. Additionally, one drop of apraclonidine or brimonidone should be instilled one hour prior to cataract surgery for all patients.
- b. Anesthetize the eye using general, retrolubar, peribulbar or topical anesthesia per standard operating procedures.
- c. Dilate the pupil using tropicamide and phenylephrine, and apply a topical non-steroidal anti-inflammatory medication.
- d. Remove the cataract using phacoemulsification techniques.
- e. Insert an intraocular lens.
- f. If the angle needs to be deepened after cataract surgery, an intracameral miotic should be instilled.

2. Stent Implantation

- a. The peel pouch containing the iStent Trabecular Micro-Bypass Stent System should be opened onto the sterile field. Caution: Do not use the device if the Tyvek® lid has been opened or the packaging appears damaged. In such cases, the sterility of the device may be compromised.
- b. The stent implantations are designed for nasal placement; therefore, surgery is recommended to be performed from the temporal side of the head.
- c. Ensure that the orientation of the stent on the inserter is appropriate for the desired nasal implantation and that the preferred hand position is on the inserter.
- d. Inspect angle with a gonioprism to ensure that a good view is available at the nasal implant location.
- e. Place a gonioscope on the cornea and reposition the surgical microscope as needed to visualize the trabecular meshwork, through the gonioprism, on the nasal side of the eye. Focus on the landmarks in the angle of the eye (**Figure 1**). Look up from the iris root to find the scleral spur (white line). Then look for Schwalbe’s line (white line) down from the cornea. The trabecular meshwork (red/brown line) is between the sclera spur and Schwalbe’s line. Schlemm’s canal is behind the trabecular meshwork.


Figure 1

3. Insertion of stent

- a. Inject viscoelastic into the anterior chamber to assist with chamber maintenance.
- b. Insert the stent (which is attached to the inserter tip) through the temporal incision that was used to extract the cataract and insert the intraocular lens.
- c. Traverse the anterior chamber with the inserter and position the inserter tip at approximately the pupillary margin. Place the gonioprism into the desired position.
- d. Locate the trabecular meshwork, and look for bifurcated areas based on asymmetric areas of pigmentation, and select an implant location below the horizontal midline of the meshwork and adjacent to any pigmented areas (which could represent collector channels).
- e. Gently slide the leading edge of the stent through the trabecular meshwork and into Schlemm’s canal at the nasal position (3 to 4:00 o’clock position for the right eye; 8 to 9:00 o’clock position for the left eye), with the tip of the implant directed inferiorly. A slight lifting motion may be required for engagement. If there is difficulty with insertion at the desired location, try inserting about 0.5 clock hour inferior. Continue to move inferior as needed for subsequent attempts. Backing up may prevent the tip of the device from penetrating tissue due to the circular geometry of the edge.
- f. Release the stent by pushing the button on the inserter. Once the stent is in Schlemm’s canal, gently tap the stent with the inserter to drop the heel of the stent completely into Schlemm’s canal. The heel of the stent will not drop into Schlemm’s canal without a gentle push.
- g. Ensure the stent is fully released and within Schlemm’s canal, and withdraw the inserter.

4. Ensuring proper placement of stent (intraoperative confirmation)

At the end of the procedure, the following should be performed:

- a. Flush the anterior chamber (AC) of any refluxed blood. Note: minimal blood reflux is a normal physiological response to placement of the stent.
- b. High-magnification examination to confirm that base of the implant is parallel with the circumferential axis of Schlemm’s canal.
- c. Nudge the inlet gently in an anterior/posterior direction to confirm that rails on the base of the implant are located on the back wall of Schlemm’s canal, and that the inlet axis is parallel with the iris plane.
- d. Nudge the inlet toward the tip of the implant in the direction of the circumferential axis of Schlemm’s canal to ensure that the inlet end of the base is well seated and fully through the trabecular meshwork.
- e. Observe with a gonioprism to verify that the inlet is radial in direction relative to the anterior chamber.
- f. Irrigate the anterior chamber with balanced salt solution (BSS) through the corneal wound to remove all viscoelastic. Press down on the posterior edge of the incision as needed to facilitate complete removal of viscoelastic.
- g. Inflate the anterior chamber with saline solution as needed to achieve physiologic pressure.
- h. Ensure that the corneal incision is sealed, and place 10-0 nylon suture *only if needed*.
- i. Prophylactic ophthalmic medications should not be administered at the end of the procedure to mitigate IOP spikes. Because a postoperative IOP check is recommended at 3 to 7 hours postoperative, the use of prophylactic medications should be implemented after the IOP is measured during the initial postoperative period. Medications can then be added if an IOP spike is observed (an IOP spike is considered to be ≥ 10 mmHg from the baseline value).

5. Postoperative Instructions

- a. At the physician’s discretion, patients may be instructed to wear an eye pad for 4 hours if retrolubar injection was done, then wear eye shield for the entire first day, and then for 1 week (at night only).
- b. Patients should be directed to administer a topical fluoroquinolone antibiotic for one week postoperatively. In addition, the patient should be instructed to administer prednisolone acetate 1% ophthalmic suspension 6 drops per day for week 1, 5 drops per day for week 2, 4 drops per day for week 3, and 3 drops per day for week 4. At the 1-month follow-up visit, the physician should determine if the patient will continue with drug tapering by reducing one drop per week, or if the medication should be discontinued.
- c. If pressure increases postoperatively requiring medical therapy, medications should be added at the physician’s discretion based upon IOP, optic nerve and visual field status. A general guidance is that medications should be considered when the unmedicated IOP reaches > 21 mmHg.
- d. If a pressure spike is observed, wound “burping” may be performed at the physician’s discretion.

ADVERSE EVENT REPORTING

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as device related and that were not previously expected in nature, severity, or incidence must be reported to Glaukos Corporation at:

Phone Number: +1.949.367.9600
Fax Number: +1.949.297.4540

HOW SUPPLIED

The stent is attached to the tip of a single-use inserter, and the system is provided sterile and nonpyrogenic in a blister tray. Each stent system is individually serialized, and the serial number is provided on the tray lid and unit carton. The device has been sterilized by gamma radiation.

STORAGE REQUIREMENTS

The device should be stored at room temperature in the range of 15-30°C.

EXPIRATION DATE

The expiration date on the device package (tray lid) is the sterility expiration date. In addition, there is a sterility expiration date that is clearly indicated on the outside of the unit carton. Sterility is assured if the tray seal is not punctured or damaged until the expiration date. This device should not be used past the indicated sterility expiration date.

RETURN GOODS POLICY

Please contact Glaukos Corporation.

CLINICAL TRIAL RESULTS

Description of the Randomized Clinical Trial

The study was a prospective, randomized, controlled, open-label, multicenter trial

conducted at 27 sites throughout the U.S. A total of 240 patients meeting the study eligibility criteria were randomized in a 1:1 fashion to undergo either implantation of the iStent in conjunction with cataract surgery (treatment group) or cataract surgery without implantation of the iStent (control group). Subjects in this randomized population were treated from April 13, 2005 through June 28, 2007. All subjects were followed for a period of 2 years. A total of 117 subjects were enrolled in the treatment group, and 123 subjects were enrolled in the control group. To obtain additional safety information, the study also included a separate non-randomized arm of an additional 46 patients who underwent iStent implantation during cataract surgery. All study subjects were diagnosed with mild to moderate open-angle glaucoma (OAG) in the study eye (including pseudoexfoliative and pigmentary glaucoma) based upon the protocol inclusion and exclusion criteria.

Of the 117 subjects in the treatment group, over two-thirds (69%) were 70 years of age or older at the time of surgery, with a mean age of 74 years. Most subjects (61%) were female, and the majority of subjects (71%) were Caucasian. There were equal proportions of right eyes and left eyes. Similar demographic characteristics were seen in the 123 control subjects, where almost two-thirds (65%) of these subjects were 70 years of age or older at the time of surgery, with a mean age of 73 years. Over half of these subjects (58%) were female. The majority of subjects were Caucasian (72%), and there were equal proportions of right eyes and left eyes.

The primary effectiveness outcome was defined as IOP ≤ 21 mmHg without use of antiglaucoma medication at 12 months. The proportion of subjects with this outcome was compared between the two study groups. The secondary effectiveness outcome was defined as IOP reduction from baseline of ≥ 20% without use of anti-glaucoma medication at 12 months. The proportion of subjects with this outcome was compared between the two study groups.

Efficacy Results – Randomized Trial

As can be seen in Table 1, 73% of subjects in the treatment group (combined cataract and iStent implantation) met the primary endpoint of IOP ≤ 21 mm Hg with no medications at 12 months. In comparison, only 50% of subjects in the control group (cataract surgery only) met the primary endpoint. This treatment difference of 23% in favor of the iStent group on the primary endpoint at 12 months was statistically (p≤ .001) and clinically significant. The secondary endpoint in the GC-003 pivotal trial was the proportion of patients with IOP reduced ≥ 20% from baseline without medications at 12 months. In the iStent treatment group, 67% of subjects implanted met this endpoint compared to 47% in the cataract control group (Table 2). This treatment difference of 20% was also statistically (p=.001) significant.

TABLE 1
IOP ≤21 mmHG WITHOUT ANTI-GLAUCOMA MEDICATIONS AT 12 MONTHS
(EXCLUDING OBSERVED OUTCOMES AFTER SECONDARY PROCEDURES)
(ITT POPULATION)

Analysis Population and Imputation Method	Cataract Surgery with iStent		Cataract Surgery Only	P-value (one-sided Z-test)
	n/N (%) (90% CI)		n/N (%) (90% CI)	
	Difference (90% CI)			
ITT Using LOCF	85/117 (73%) (66%,79%)		61/123 (50%) (42%,57%)	<.001
	23% (13%, 33%)			

Secondary procedures include stent reposition, stent explant/replacement, 2 stents implanted, deep sclerectomy, replaced IOL, trabeculoplasty, pars plana vitrectomy, and iridoplasty for stent obstruction. Observed IOP after secondary procedures was treated as missing. Subjects with anti-glaucoma medications were treated as a non-responder regardless of receiving with secondary surgical intervention.

TABLE 2
IOP REDUCTION ≥ 20% WITHOUT ANTI-GLAUCOMA MEDICATIONS AT 12 MONTHS
(EXCLUDING OBSERVED OUTCOMES AFTER SECONDARY PROCEDURES)
(ITT POPULATION)

Analysis Population and Imputation Method	Cataract Surgery with iStent n/N (%) (90% CI)	Cataract Surgery Only n/N (%) (90% CI)	P-value (one-sided Z-test)
	Difference (90% CI)		
	ITT Using LOCF	78/117 (67%) (59%,74%) 20% (9%, 30%)	
			0.001

Secondary procedures include stent reposition, stent explant/replacement, 2 stents implanted, deep sclerectomy, replaced IOL, trabeculoplasty, pars plana vitrectomy, and iridoplasty for stent obstruction. Observed IOP after secondary procedures was treated as missing. Subjects with anti-glaucoma medications were treated as a non-responder regardless of receiving with secondary surgical intervention.

Safety Results – Adverse Events

A summary of ocular adverse events reported in the safety population during the randomized clinical trial is presented in Table 3. The Safety population includes all subjects who underwent cataract surgery, with or without the iStent procedure, according to the actual treatment each subject received. This population consists of 233 subjects (111 treatment and 122 control).

TABLE 3
OCULAR ADVERSE EVENTS
SAFETY POPULATION

Adverse Events	Cataract Surgery with iStent N = 111 n (%)	Cataract Surgery Only N = 122 n (%)
Secondary surgical intervention (including intervention for elevated IOP)	37 (33%)	39 (32%)
Increase in IOP: ≥ 10 at any time ≥ at or after 1 month postoperative	34 (31%) 0 (0%)	38 (31%) 2 (2%)
Any BCVA loss of at least 1 line at or after the three month visit	7 (6%)	9 (7%)
Elevated IOP requiring treatment with oral or intravenous medications	3 (3%)	2 (2%)
Epiretinal membrane	2 (2%)	1 (1%)
Age related macular degeneration	2 (2%)	0 (0%)
Iris atrophy	2 (2%)	0 (0%)
Stent obstruction by iris, vitreous, fibrous overgrowth, fibrin, blood, etc.	2 (2%)	0 (0%)
Iritis	1 (1%)	5 (4%)
Allergic conjunctivitis	1 (1%)	2 (2%)
Corneal edema	1 (1%)	2 (2%)
Disc hemorrhage	1 (1%)	2 (2%)
Dry eye	1 (1%)	2 (2%)
Cystoid Macular Edema	1 (1%)	1 (1%)
Macular edema	1 (1%)	1 (1%)
Chronic pain due to medication (Xalatan)	1 (1%)	0 (0%)
Corneal abrasion	1 (1%)	0 (0%)
Endopigment	1 (1%)	0 (0%)
Eye splash injury	1 (1%)	0 (0%)
Metallic particle on iris	1 (1%)	0 (0%)
Possible bacterial conjunctivitis	1 (1%)	0 (0%)
Subconjunctival hemorrhage	1 (1%)	0 (0%)
Uveitis	1 (1%)	0 (0%)
Posterior capsular opacification	0 (0%)	3 (2%)
Anterior chamber +1 cells (at 1 month) requiring treatment	0 (0%)	2 (2%)
Blurry vision	0 (0%)	2 (2%)
Foreign body sensation	0 (0%)	2 (2%)
Blepharoconjunctivitis	0 (0%)	1 (1%)
Carotid artery disease with loss of BCVA	0 (0%)	1 (1%)
Choroidal tubercle	0 (0%)	1 (1%)
Conjunctival injection	0 (0%)	1 (1%)
Conjunctivitis	0 (0%)	1 (1%)
Endophthalmitis	0 (0%)	1 (1%)
Episcleritis	0 (0%)	1 (1%)
Intermittent tearing	0 (0%)	1 (1%)
Keratitis sicca	0 (0%)	1 (1%)
Lesion on eyelid	0 (0%)	1 (1%)
Macular hole	0 (0%)	1 (1%)
Macular traction	0 (0%)	1 (1%)
Posterior capsule bag rupture	0 (0%)	1 (1%)
Postoperative refractive error	0 (0%)	1 (1%)
Proliferative diabetic retinopathy	0 (0%)	1 (1%)
Segmental loss of neuroretinal rim	0 (0%)	1 (1%)
Superficial punctate keratitis	0 (0%)	1 (1%)

Conclusions

The iStent® Trabecular Micro-Bypass Stent Model GTS100R/L was found to be safe and effective in conjunction with cataract surgery for the reduction of IOP in subjects with mild to moderate open-angle glaucoma who were being treated with ocular hypotensive medication.

MRI SAFETY INFORMATION:

Static Magnetic field of 3 Tesla or less

Non-clinical testing has demonstrated that the iStent Trabecular Micro-Bypass Stent (Models GTS100R and GTS100L) is MR Conditional.

A patient with this device can be safely scanned in an MR system meeting the following conditions:
· Static magnetic field of 3T or less
· Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m)
· Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)

Under the scan conditions defined above, the iStent Trabecular Micro-Bypass Stent (Models GTS100R and GTS100L) is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends less than 15 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Static magnetic field of 7-Tesla, ONLY:

Non-clinical testing demonstrated that the iStent Trabecular Micro-Bypass Stent (Models GTS100R and GTS100L) is MR Conditional.

A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:
· Static magnetic field of 7-Tesla, ONLY
· Maximum spatial gradient magnetic field of 10,000-gauss/cm (extrapolated) or less
· Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
· First Level Controlled Operating Mode of operation for the MR system
· Use of a transmit/receive RF head coil, ONLY

MRI-Related Heating

In non-clinical testing, the device produced a temperature rise of 0.4°C during MRI performed for 15-minutes of scanning (i.e., per pulse sequence) in a 7-Tesla/298-MHz MR system (Philips Acheiva, Philips Healthcare, Cleveland, OH) using at transmit/receive RF head coil.

MR system reported, whole body averaged: SAR < 1-W/kg

Artifact Information

In non-clinical testing, the image artifact caused by the device extends approximately 23-mm from this device when imaged using a gradient echo pulse sequence and a 7-Tesla MR system.

LABELING

The following symbols are used on the device packaging.

<i>Symbol</i>	<i>Description</i>
	Catalogue/Model Number
	Serial Number (for the stent)
	Lot Number
	Do not re-use
	Use By (year-month-day)
	Do not use if package is damaged
	Sterilized by Gamma Irradiation
	Temperature Storage Requirements
	CE Marking
	MR Conditional
	For prescription use only
	Consult Instructions For Use.
	Manufacturer
	Authorized European Representative

CAUTION

U.S. Federal law restricts this device to sale by, or on the order of, a physician.

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