
PRESERFLO™ MicroShunt **Instructions for Use (en)**

INTRODUCTION

The PRESERFLO MicroShunt product consists of the PRESERFLO MicroShunt and a 3-20mm stainless steel Scleral Marker. The MicroShunt is an implantable glaucoma drainage device made of an extremely flexible SIBS [poly (Styrene-block-IsoButylene-block-Styrene)] polymer with a tube of 350 µm outer diameter and a lumen of 70 µm. It has triangular fins that prevent migration of the tube into the AC. The device is designed to be implanted under the subconjunctival/Tenon space.


1. Make sure to check the expiration date on packaging

INTENDED USE/ PURPOSE

The PRESERFLO MicroShunt employs a tube to create a conduit for the flow of aqueous humor from the anterior chamber of the eye to a bleb formed under the conjunctiva and Tenon's capsule; the front/ proximal end of the tube extends into the anterior chamber while the back/ distal end terminates in the bleb. The PRESERFLO MicroShunt reduces IOP by physically shunting aqueous from the high pressure anterior chamber to the lower pressure bleb.

PRESERFLO™ MicroShunt CONTENTS (QTY.)

- Case Report Labels (12)
- MicroShunt (1)

- Patient Card (1)
- 3mm Scleral Marker (1)

RECOMMENDED ACCESSORIES (QTY.)

- The PRESERFLO™ MicroShunt recommended accessories for surgical procedure depending on the configuration are:
- Marker Pen – Gentian Violet (1)
 - Anterior Chamber Cannula 23G 8mm bend (1)
 - MANI Ophthalmic Knife Slit-Angled 1.0mm Knife (1) or, Ophthalmic Knife Double Step-Angled 1.0mm Knife (1)
 - Sponges (3) and
 - Sclera Track Needle 25g x 5/8 (25G Needle) (1)

INTENDED USERS

The intended users of the PRESERFLO MicroShunt are ophthalmologists/ophthalmic surgeons specializing in the treatment of glaucoma (including surgeons specializing in anterior segment and cataract surgery) who have been trained to use the device.

INDICATIONS FOR USE

The MicroShunt is intended for reduction of intraocular pressure in eyes of patients with primary open angle glaucoma where IOP remains uncontrollable while on maximum tolerated medical therapy and/or where glaucoma progression warrants surgery.

INTENDED PATIENT/ TARGET GROUPS

The target patient population is adult patients with primary open-angle glaucoma where IOP remains uncontrolled while on maximum tolerated medical therapy and/or where glaucoma progression warrants surgery.

CONTRAINDICATIONS

The implantation of the MicroShunt is contraindicated under the following circumstances and conditions: ANGLE CLOSURE GLAUCOMA; PRESENCE OF CONJUNCTIVAL SCARRING, PREVIOUS INCISIONAL OPHTHALMIC SURGERY INVOLVING THE CONJUNCTIVA OR OTHER CONJUNCTIVAL PATHOLOGIES (E.G., THIN CONJUNCTIVA, PTERYGIUM) IN THE TARGET QUADRANT; ACTIVE IRIS NEOVASCULARIZATION; ACTIVE INFLAMMATION (E.G., BLEPHARITIS, CONJUNCTIVITIS, SCLERITIS, KERATITIS, UVEITIS); VITREOUS IN THE ANTERIOR CHAMBER; PRESENCE OF AN ANTERIOR CHAMBER INTRAOCULAR LENS (ACIOL); INTRAOCULAR SILICONE OIL.

WARNINGS

- Rx only: This device is restricted to sale by, or on the order of, a physician.
- For one-time use only. Do not reuse or re-sterilize. Reuse, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in serious patient injury, illness, blindness or death. Reuse, or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, blindness, or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.
- Long term effects of Mitomycin C (MMC) with the use of this device have not been evaluated. Necessary precautions and interventions on the use of MMC are highly recommended.
- The MicroShunt should not be subjected to direct contact with petrolatum-based (i.e., petroleum jelly) materials (e.g., ointments, dispersions, etc.).

- The effects of cutting or modifying the MicroShunt have not been evaluated.
- Viscoelastics have not been tested with this device. However, in an emergency when all other therapies have failed, the use of hydroxypropyl methylcellulose (HPMC) may be an option. Use of HPMC should be a last resort to correct a flat chamber with the MicroShunt and may risk loss of flow through the device for one or more weeks after use necessitating close or more frequent observation of IOP.

PRECAUTIONS

- The safety and effectiveness of the MicroShunt has not been established in patients with the following circumstances or conditions:
 - o Eyes with chronic inflammation
 - o Eyes with congenital and infantile glaucoma
 - o Eyes with neovascular glaucoma
 - o Eyes with uveitic glaucoma
 - o Eyes with pseudoexfoliative or pigmentary glaucoma
 - o Eyes with other secondary open angle glaucomas
 - o Eyes that have undergone previous incisional glaucoma surgery or cilioablativ procedures.
 - o Implantation with concomitant cataract surgery with IOL implantation
- Avoid use of toothed forceps to handle device. McPherson type forceps are recommended.
- The patient's IOP should be monitored postoperatively. If IOP is not adequately maintained after surgery, appropriate additional therapy to maintain IOP should be considered.
- The safety and effectiveness of the use of more than a single PRESERFLO™ MicroShunt has not been established.
- If the MicroShunt appears deformed, folded and/or distorted, do not use.
- The creation of 2 or more scleral tunnels in an attempt to implant the MicroShunt may cause leakage of aqueous humor and increase the risk of hypotony, if the additional scleral tunnels did not self-seal or sutured.
- When managing patients with a shallow or flat anterior chamber, while not a likely occurrence, there is a possibility that pharmacologic dilation of the pupil in these patients may cause the iris to obstruct the MicroShunt. For these patients, more frequent monitoring of IOP following dilation is recommended.

POTENTIAL COMPLICATIONS/ADVERSE EVENTS

The complications during and after surgery may include:

Glaucoma progression not controlled, increase in cup-to-disc ratio (C/D), anesthesia related complications, difficulty in inserting the MicroShunt, failure to implant the device, device malfunction, device repositioning, extended surgical procedure, tube migration out of anterior chamber, flat anterior chamber, shallow anterior chamber, excessive bleeding in anterior chamber or incision site or eye, MicroShunt touches cornea or iris, intraocular pressure too high or low, strabismus, viscoelastic used in anterior chamber, choroidal effusion or hemorrhage, retinal complications (retinal detachment, proliferative retinopathy), hyphema, hypotony or hypotony maculopathy, phthisis bulbi, endophthalmitis, tube erosion through conjunctiva, tube obstruction partial or complete (block by iris or vitreous or fibrin or debris), iridodialysis, iris incarceration, uveitis, iritis, diplopia, aqueous misdirection, corneal complications (abrasion, edema, ulceration, infection, decompensation, bullous keratopathy, endothelial cell loss, Descemet striae, keratitis), partial or complete vision loss, blurry vision, bleb leak, blebitis, cystic bleb, bleb failure, pupillary block, ptosis, macular edema, prolonged inflammation, use of glaucoma medications, ocular pain, conjunctival complications (buttonhole, dehiscence, dissection, hemorrhage, hyperemia, scar, tear, ulcer), iris adhesions/synechiae, iris abrasion, cataract development or progression, posterior capsule opacity, explantation of the MicroShunt, encapsulation reaction, foreign body sensation, fibrin in anterior chamber, visual field damage, unplanned glaucoma-related surgical re-intervention, optic disc hemorrhage, globe perforation, headache, vitreous hemorrhage, wound leak, and suture related complications.

STERILIZATION

All components of the PRESERFLO™ MicroShunt are sterile upon receipt. The components of the PRESERFLO™ MicroShunt are sterilized by ETO sterilization per ISO 11135 guidelines. Contents will remain sterile unless the package is opened or damaged. Should the PRESERFLO™ MicroShunt be inadvertently rendered non-sterile, or should the package be opened and/or damaged, discard the device. **DO NOT RESTERILIZE.**

STORAGE

Store at controlled room temperature (20±5°C) and avoid excessive heat, excessive humidity, direct sunlight, and water.

INSTRUCTIONS FOR USE

User must be properly trained on the use of the device prior to use.

Surgery:

- Anesthesia:** The type of anesthesia to be administered during the procedure is at the surgeon's discretion.
- Preparation of PRESERFLO™ MicroShunt:** a) Carefully examine the pouch containing the PRESERFLO™ MicroShunt. If the pouch is damaged, do not use the device. b) Remove PRESERFLO™ MicroShunt from sterile packaging onto sterile field. c) Wet the MicroShunt using a solution of Balanced Salt Solution (BSS).
- Planning of Implant Position:** Insertion of the MicroShunt is at the surgeon's discretion, either in the superonasal or superotemporal quadrant of the eye, avoiding the rectus muscles.
- Corneal Traction and Conjunctival/Tenon's Flap:** Perform corneal traction with suture at surgeon's discretion. In the planned surgical area, make a 6-8mm peritomy incision, then dissect a subconjunctival/subTenon pouch (fornix based) over a circumference of 90°-120° from near the limbus to at least 8mm deep. Assure there is no active bleeding in the area of the planned MicroShunt entry. Avoid excessive use of cautery, which could increase the risk of scarring.
- Application of MMC:** At the surgeon's discretion, apply sponges saturated with MMC on the sclera under the subconjunctival/subTenon pouch, making sure they do not come into contact with the conjunctival wound edge. Remove the sponges and rinse with generous amounts of BSS (>20mL).
- Preparation of Anterior Chamber Entrance Site:** Ink the 3-20mm Marker (Scleral Marker) points with the Marker Pen - Gentian Violet and mark a point 3mm from the limbus.
 - If using the **Double Step-Angled 1.0mm Knife:** At the marked point, create a scleral tunnel into the anterior chamber (AC) entering just above the iris plane (level of trabecular meshwork). As a guide externally, enter the AC just below the middle border of the surgical limbus (posterior border of blue/gray zone, which should overlie internally with Schwalbe's line), if present. The scleral tunnel is complete when a scleral pocket that is 2mm long is created at the distal end of the scleral tunnel, where the fin portion of the MicroShunt is tucked.
 - If using the **Slit Angled 1.0mm Knife:** At the marked point, create a shallow scleral pocket that is 2mm long, where the fin portion of the MicroShunt is tucked. Insert a 25G needle into the scleral pocket and create a scleral tunnel into the AC entering just above the iris plane (level of trabecular meshwork). As a guide externally, enter the AC just below the middle border of the surgical limbus (posterior border of blue/gray zone, which should overlie internally with Schwalbe's line), if present.
- Insertion of MicroShunt into Anterior Chamber: Using a non-toothed forceps,** position the MicroShunt tube with the bevel facing towards the cornea close to the scleral pocket/scleral tunnel while keeping the fin face parallel to the scleral surface. Gently insert the MicroShunt into the scleral pocket/scleral tunnel until the fin is tucked under the scleral pocket. Use of excessive force will result in temporary bending of the MicroShunt and increase the difficulty of entrance. If not accessible, another scleral pocket/scleral tunnel can be made approximately 1mm to either side of the original tract. If the MicroShunt is damaged during the procedure, use a new MicroShunt and use a new scleral pocket/scleral tunnel. After successful insertion of the MicroShunt, check that the proximal end of the tube in the AC is not in contact with the cornea or iris, and that there is flow of aqueous from the distal end. Then tuck the end of the MicroShunt tube under the subconjunctival/subTenon's flap on the scleral surface, making sure it is straight and not caught in Tenon's capsule. If flow is not visible, consider one or more of the following actions: applying gentle pressure on the eye (sclera, limbus, cornea), filling the anterior chamber with BSS through a paracentesis, or using a 23G thin wall 8mm bent cannula inject BSS into the distal end of MicroShunt to prime the tube. If flow is not established, reposition, or replace the device.
- Verify position of the MicroShunt in the anterior chamber angle:** It is recommended that intraoperative gonioscopy be performed to assess the position of the proximal tip of the MicroShunt in the anterior chamber to ensure that there is sufficient clearance from the cornea or iris.
- Conjunctival Closure:** Reposition the Tenon's capsule and conjunctiva to the limbus and suture using sutures with a long history of successful use by glaucoma specialists (i.e., nylon, vicryl). A moistened fluorescein strip is used to check for a conjunctival leak, which may be closed with additional sutures. Verify the following: presence of the proximal end of the MicroShunt in the anterior chamber; distal end of the MicroShunt is straight and not caught in Tenon's capsule. Pull corneal traction suture out if it was used. Apply antibiotic and steroidal medication postoperatively, as appropriate.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy. Follow the MMC instructions for use for proper disposal of the MMC.

Post operatively:

- Monitor the intraocular pressure at each follow-up to determine if the MicroShunt is patent.
- If the MicroShunt is repositioned, removed and/or replaced with another type of device based on the device performance, the conjunctiva should be cut at the limbus in a similar manner to the original procedure. The device should be exposed, and the repositioning or removal conducted. If the device is removed, verify that there is no aqueous leakage from the track. If there is, use a suture the track to seal it. Conjunctival closure should follow the original procedure.

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Santen UK Ltd
UK Responsible Person Santen UK Ltd Salisbury Hall, St Albans, AL2 1BU, United Kingdom