Micro-Bypass System		P-00588
INSTRUCTIONS FOR USE		Rev 4/03/26
DEVIce	DESCRIPTION	
Catalogue #	Description	
153	Glaukos® iStent Infinite® System Model IS3 with 3 (3) G2-W stearalkonium heparin coated Trabecular Micro-Bypass Stents preloaded in one injector	

The iStent infinite® System contains (3) intraocular stents that are manufactured from implant grade titanium (Ti6Al4V ELI) and are coated with stearalkonium heparin (note: the heparin is from a porcine source). The stent has a single piece design, is approximately 360 microns in height and 360 microns in diameter. The injector has been designed by Glaukos Corporation to implant three Model G2-W stents through a single incision into Schlemm's canal.

The implant is designed to stent open a passage through the trabecular meshwork to allow for an increase in the facility of outflow and a subsequent reduction in intraocular pressure.

Each iStent Infinite System is provided sterile in a blister tray (only the inner contents of the blister tray are sterile). Three stents are already assembled onto the single-use injector. The blister tray lid is labeled with the required product identification information. A lot number and expiration date are provided on the labeling. The iStent infinite System is sterilized by gamma irradiation.

The expiration date is a sterility expiration date (tryd) is the sterility expiration date. In addition, there is a device expiration date that is clearly indicated on the outside of the unit carton. Sterility is assured if the tray seal is not broken, punctured or damaged until the expiration date. This device should not be used past the indicated sterility expiration date.

INDICATIONS FOR USE

The iStent infinite, is intended to reduce intraocular pressure safely and effectively in adult patients diagnosed with primary open-angle glaucoma, pseudo-exfoliative glaucoma or pigmentary glaucoma.

The device is safe and effective when implanted in combination with or without cataract surgery, or in combination with glaucoma pressure reduction and/or would benefit from glaucoma medication reduction.

The device may also be implanted in patients who continue to have elevated intraocular pressure despite prior treatment with glaucoma medications and/or conventional glaucoma surgery.

CONTRAINDICATIONS

The iStent infinite System is contraindicated under the following circumstances or conditions:

- In eyes with primary angle closure glaucoma, or secondary angle-closure glaucoma, including neovascular glaucoma, because the device would not be expected to work in such situations
- In patients with retinal or vitreous disease, Sturge-Weber Syndrome or any other type of condition that may cause elevated episcleral venous pressure

INSTRUCTIONS FOR USE

- Remove the Stent Delivery Button Lock from the injector.
- Make a corneal incision to allow entry of the injector into the anterior chamber. Recommended incision location is the temporal peripheral cornea for either eye. Ophthalmic viscosilastic (cohesive) should be used to form the anterior chamber, as necessary. Deepen the anterior chamber by injecting with viscosilastic as needed, being careful not to overfill.

The iStent infinite injector insertion steps are as follows: With the gonioscope removed from the cornea, insert the injector introducer Tip through the clear corneal incision into the anterior chamber. Once the introducer Tip has been inserted through the corneal incision, the Introducer Tip auto-retracts exposing the Insertion Tube and Trocar tip. Advance the Insertion Tube to the pupillary margin toward the targeted trabecular meshwork tissue (i.e., the ob/ interno approach). Take care to avoid contact with the IOL, iris, or cornea.

Place the gonioscope on the cornea and position the patient and surgical microscope as needed to visualize the trabecular meshwork through the gonioscope on the nasal side of the eye.

Advance the Insertion Tube containing the trocar towards the TM (just above the scleral spur) and penetrate the trabecular meshwork at the center of the TM. The trocar is used to not only penetrate the TM, but will remain in the tissue to act as an axial guide for the stent as the stent traverses over the trocar through to Schlemm's canal.

Gently hold the Insertion Tube against the TM and apply appropriate pressure to slightly indent or "dimple" the tissue (those should stretch just enough to form a "V" when pressing on the TM).

Slowly squeeze and hold down the Stent Delivery Button to automatically inject the stent head through the TM and into Schlemm's canal. Look through the window in the insertion tube to verify that the stent has been implanted properly.

Withdraw the injector straight back from the stent, and verify that the stent is well-positioned and secured in the TM.

Actuate the Singulator to prepare the next stent for implantation. Listen for two audible clicks, and verify that the next stent is visible just above the scleral spur.

Repeat the previous step for the third stent, placing the third stent approximately 2 clock hours away from either of the first two stents.

Remove the injector from the eye.

Irrigate the anterior chamber with balanced salt solution (BSS) through the corneal wound to remove all viscosilastic. Press down on the posterior edge of the cornea as needed to facilitate complete removal of viscosilastic.

Inflate the anterior chamber with saline solution as needed to achieve normal physiologic pressure.

Ensure that the corneal incision is sealed.

Relative contraindications

- Prep the patient as one would for stent implantation surgery.
- Re-open the eye at the preferred location in order to reach the stent. A clear corneal incision measuring approximately 1.5mm in length is recommended.
- Use viscosilastic to inflate the anterior chamber to create access to the stent's location, move the stent away from a delicate structure if loose, and/or protect intraocular tissues.
- Use a gonioscope if needed to visualize the location of the stent in the anterior chamber.
- Use a forceps device through the corneal incision and grasp the stent in a convenient and secure manner before removing the stent from the anterior chamber.
- Irrigate the anterior chamber with balanced salt solution (BSS) through the corneal wound to remove all viscosilastic. Press down on the posterior edge of the cornea as needed to facilitate complete removal of viscosilastic. Repeat as needed until all viscosilastic has been removed.

Static magnetic field of 3-Tesla or less:

Non-clinical testing has demonstrated that the iStent infinite Trabecular Micro-Bypass System, Model IS3 is MR Conditional.

MR

A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field: 3T or less
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg (First Level Controlled Operating Mode)

Under the scan conditions defined above, the iStent infinite Trabecular Micro-Bypass System, Model IS3 is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends less than 15 mm from the device when imaged with a gradient echo pulse sequence and a 3T MRI system.

Static magnetic field of 7-Tesla, ONLY:

Non-clinical testing demonstrated that the iStent infinite Trabecular Micro-Bypass System, Model IS3 is MR Conditional.

MR

A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- Static magnetic field: 7T, ONLY
- Maximum spatial gradient magnetic field of 10,000-gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
- First Level Controlled Operating Mode of operation for the MR system
- Use of a transmit/receive RF head coil, ONLY

MR-Related Heating

In non-clinical testing, the device produced a temperature rise of 0.4°C during MRI performed for 15-minutes of scanning (i.e., per pulse sequence) at a 7-Tesla MRI system (Philips Achieva, Philips Healthcare, Cleveland, OH) using a transmit/receive RF head coil. MR system reported, whole body averaged: SAR < 1-W/kg.

Artifact Information

In non-clinical testing, the image artifact caused by the device extends approximately 23-mm from this device when imaged using a gradient echo pulse sequence and a 7T MRI system.

LABELING

The following symbols are used on the device packaging.

Symbol **Definition**

REF Catalogue/Model Number

SN Serial Number (for the stent)

LOT Lot Number

YYMMDD Do not reuse

YYMMDD Use By (year-month-day)

Do not use if package is damaged

STERILE | R Sterilized by Gamma Irradiation

TEMPERATURE STORAGE REQUIREMENTS Temperature Storage Requirements

MR MR Conditional

Rx Only For Prescription use only

i Consult Instructions For Use

Manufacturer Manufacturer

QTY Quantity

MD Medical Device

Single sterile barrier system with protective packaging outside

Glaukos® and iStent infinite® are registered trademarks of Glaukos Corporation.

Typex® is a registered trademark of DuPont USA.

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Gaukos® dan iStent infinite® adalah merek dagang terdaftar milik Glaukos Corporation.